

Richmond, April 26th 2018

Product Updates

Datrend Systems Inc. is pleased to announce the following updates to the vPad range of products and solutions:

Printing of PM Labels

The latest release of vPad Apps features wireless printing of PM labels for inspected equipment. Compatible with a range of Brother printers, communication can be via a WiFi network or directly via Bluetooth.

Improved Accuracy of vPad-RF

V1.32 of the vPad-RF App was released on the 27th of March 2018. vPad-RF will now provide multi-pulse analysis, which averages measurements from multiple RF pulses and reduces errors caused by analyzing pulses in isolation. This update further improves the accuracy of this unique product.

Release of vPad-Medimizer

The vPad-Medimizer App communicates with MediMizer (version 14+) via the internet to download PM schedules, work orders, and equipment information, and upload completed Test Records to the CMMS, enabling work orders to be closed automatically. Used in conjunction with vPad-Check, this methodology eliminates the need for paper records thereby improving throughput, accuracy and traceability.

Masimo Rainbow

The vPad-A1 App has been updated to enable routine preventive maintenance of patient monitoring systems which incorporate Rainbow technology. A White Paper has been produced detailing how the technology works and instruction on how to test rainbow monitors using the vPad-A1 All-In-One Vital Signs Tester. This white paper is available on our website at www.datrend.com/support/white-papers/

For more information on any of these updates visit the company's website www.datrend.com or contact them directly at sales@datrend.com

About Datrend Systems Inc.

Datrend Systems Inc. is a Canadian based biomedical engineering company focused on developing and manufacturing leading edge biomedical instruments and solutions for hospitals, clinics, third party service organizations, and medical device manufacturers. Utilizing the latest technology, we develop comprehensive Preventive Maintenance Solutions to help our customers save time and money whilst meeting stringent regulatory requirements.

Contact Information
Datrend Systems Inc.
Jon Strong
General Manager
+1-604-291-7747